

USP44 antibody



Catalog Number: 116683

Product name

USP44 antibody

Specificity

Human; other species not tested.

Antibody description

USP44 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse testis tissue. Observed molecular weight by Western-blot: 81kd

Preparation

This antibody was obtained by immunization of USP44 recombinant protein (Accession Number: NM_032147). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116683(USP44 antibody) at dilution of 1:500