

UT-B antibody



Catalog Number: 116693

Product name

UT-B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

UT-B Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 48-55 kDa

Preparation

This antibody was obtained by immunization of UT-B recombinant protein (Accession Number: NM_015865). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116693(SLC14A1 Antibody) at dilution of 1:2000