

# VAMP4 antibody

Catalog Number: 116707

## Product name

VAMP4 antibody

ELISA, IF, WB

## Specificity

Human, Mouse, Rat; other species not tested.

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Antibody description

VAMP4 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue, SH-SY5Y cells. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 18 kDa

## Validations



## Preparation

This antibody was obtained by immunization of VAMP4 recombinant protein (Accession Number: NM\_001185127). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116707(VAMP4 Antibody) at dilution of 1:600

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:116707(VAMP4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)