

VPS28 antibody

Catalog Number: 116779

Product name

VPS28 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS28 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 28-30 kDa

Preparation

This antibody was obtained by immunization of VPS28 recombinant protein (Accession Number: NM_016208). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116779(VPS28 Antibody) at dilution of 1:600



IP Result of anti-VPS28 (IP:Catalog No:116779, 4ug; Detection:Catalog No:116779 1:500) with mouse brain tissue lysate 4400ug.