

VPS37B antibody



Catalog Number: 116786

Product name

VPS37B antibody

Specificity

Human; other species not tested.

Antibody description

VPS37B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, COLO 320 cells. Observed molecular weight by Western-blot: 31 kDa

Preparation

This antibody was obtained by immunization of VPS37B recombinant protein (Accession Number: NM_024667). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116786(VPS37B Antibody) at dilution of 1:1000