

VPS53 antibody

Catalog Number: 116794

Product name

VPS53 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS53 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 95-97 kDa

Preparation

This antibody was obtained by immunization of VPS53 recombinant protein (Accession Number: BC029560). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116794(VPS53 antibody) at dilution of 1:300



IP Result of anti-VPS53 (IP:Catalog No:116794, 4ug; Detection:Catalog No:116794 1:500) with mouse liver tissue lysate 4000ug.