

VRK2 antibody

Catalog Number: 116798

Product name

VRK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VRK2 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, BxPC-3 cells, human liver tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of VRK2 recombinant protein (Accession Number: XM_005264540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116798(VRK2 antibody) at dilution of 1:500



IP Result of anti-VRK2 (IP:Catalog No:116798, 4ug; Detection:Catalog No:116798 1:500) with K-562 cells lysate 3600ug.