

# VTCN1 antibody



Catalog Number: 116806

## Product name

VTCN1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VTCN1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A549 cells, HeLa cells, rat liver tissue. Observed molecular weight by Western-blot: 55-65 kDa

## Preparation

This antibody was obtained by immunization of VTCN1 recombinant protein (Accession Number: NM\_001253849). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116806(VTCN1 Antibody) at dilution of 1:600