

# WEE1 antibody

Catalog Number: 116819

## Product name

WEE1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

WEE1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 72 kDa

## Preparation

This antibody was obtained by immunization of WEE1 recombinant protein (Accession Number: NM\_003390). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116819(WEE1 Antibody) at dilution of 1:300



IP Result of anti-WEE1 (IP:Catalog No:116819, 4ug; Detection:Catalog No:116819 1:300) with mouse brain tissue lysate 3680ug.