

WNK4 antibody



Catalog Number: 116833

Product name

WNK4 antibody

Specificity

Human; other species not tested.

Antibody description

WNK4 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, mouse liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of WNK4 recombinant protein (Accession Number: NM_032387). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116833(WNK4 Antibody) at dilution of 1:600