

# WNK4 antibody



Catalog Number: 116833

## Product name

WNK4 antibody

## Specificity

Human; other species not tested.

## Antibody description

WNK4 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, mouse liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of WNK4 recombinant protein (Accession Number: NM\_032387). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116833(WNK4 Antibody) at dilution of 1:600