

WDR19/IFT144 antibody



Catalog Number: 116858

General Information

Product name:

WDR19/IFT144 antibody

Immunogen:

Applications:

ELISA, WB, IF, IHC, IP

Specificity:

Human, Mouse, Rat; other species not tested.

Antibody description:

WDR19/IFT144 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells. Positive IHC detected in human prostate hyperplasia tissue, human prostate cancer tissue. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 145-152kd

Properties

Preparation:

This antibody was obtained by immunization of WDR19/IFT144 recombinant protein (Accession Number: BC032578). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type: Rabbit IgG

Applications

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116858(WDR19 antibody) at dilution of 1:600



IF result (ocalization near the Fleming body) of anti-WDR19 (Catalog No:116858, 1:50) with telophase hTERT-RPE1 cells (MeOH fixed) by Dr. Moshe Kim.


WDR19/IFT144 antibody

Catalog Number: 116858



IP Result of anti-WDR19 (IP:Catalog No:116858, 4ug; Detection:Catalog No:116858 1:500) with HEK-293 cells lysate 3200ug.

Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:116858 (WDR19 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:116858 (WDR19 Antibody) at dilution of 1:50

