

ZC3H11A antibody



Catalog Number: 116922

Product name

ZC3H11A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZC3H11A Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of ZC3H11A recombinant protein (Accession Number: NM_001319238). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:116922(ZC3H11A Antibody) at dilution of



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:116922(ZC3H11A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)