

# ZMYND12 antibody



Catalog Number: 116942

## Product name

ZMYND12 antibody

## Immunogen

[Human ZMYND12 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZMYND12 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 42-45 kDa

## Preparation

This antibody was obtained by immunization of ZMYND12 recombinant protein (Accession Number: NM\_032257). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116942(ZMYND12 Antibody) at dilution of 1:600