

# ZNF253 antibody



Catalog Number: 116966

## Product name

ZNF253 antibody

## Immunogen

[Human ZNF253 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZNF253 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, SMMC-7721 cells. Observed molecular weight by Western-blot: 58-65 kDa

## Preparation

This antibody was obtained by immunization of ZNF253 recombinant protein (Accession Number: NM\_021047). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116966(ZNF253 Antibody) at dilution of 1:300