

XBP1S antibody



Catalog Number: 116983

Product name

XBP1S antibody

Specificity

Human; other species not tested.

Antibody description

XBP1S Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 55 kDa, 30 kDa

Preparation

This antibody was obtained by immunization of XBP1S recombinant protein (Accession Number: NM_005080). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:116983(XBP1 Antibody) at dilution of 1:300