

# ZNF549 antibody



Catalog Number: 117003

## Product name

ZNF549 antibody

## Immunogen

[Human ZNF549 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZNF549 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, A549 cells, HeLa cells, SMMC-7721 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of ZNF549 recombinant protein (Accession Number: NM\_001199295). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:117003(ZNF549 Antibody) at dilution of 1:600



Immunofluorescent analysis of A549 cells using Catalog No:117003(ZNF549 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)