

ZNF572 antibody



Catalog Number: 117006

Product name

ZNF572 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF572 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A549 cells, HeLa cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of ZNF572 recombinant protein (Accession Number: NM_152412). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117006(ZNF572 Antibody) at dilution of 1:1000