

ZNF671 antibody



Catalog Number: 117020

Product name

ZNF671 antibody

Immunogen

[Human ZNF671 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF671 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 61-66 kDa

Preparation

This antibody was obtained by immunization of ZNF671 recombinant protein (Accession Number: BC025728). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117020(ZNF671 antibody) at dilution of 1:500