

# ZNF695 antibody



Catalog Number: 117022

---

## Product name

ZNF695 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF695 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of ZNF695 recombinant protein (Accession Number: BC023527). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117022(ZNF695 Antibody) at dilution of 1:300