

ZNF699 antibody

Catalog Number: 117023

Product name

ZNF699 antibody

Specificity

human; other species not tested.

Antibody description

ZNF699 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 74kd

Preparation

This antibody was obtained by immunization of ZNF699 recombinant protein (Accession Number: NM_198535). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117023(ZNF699 antibody) at dilution of 1:500