

ZNF271 antibody



Catalog Number: 117153

Product name

ZNF271 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF271 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of ZNF271 recombinant protein (Accession Number: BC017710). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:117153(ZNF271 antibody) at dilution of 1:1000