

ZNF354A antibody



Catalog Number: 117166

Product name

ZNF354A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNF354A Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Observed molecular weight by Western-blot: 69 kDa

Preparation

This antibody was obtained by immunization of ZNF354A recombinant protein (Accession Number: XM_047417644). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:117166(ZNF354A antibody) at dilution of 1:1000