

ZNF385D antibody



Catalog Number: 117169

Product name

ZNF385D antibody

Specificity

human, mouse; other species not tested.

Antibody description

ZNF385D Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, human testis tissue. Observed molecular weight by Western-blot: 42kd

Preparation

This antibody was obtained by immunization of ZNF385D recombinant protein (Accession Number: NM_024697). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117169(ZNF385D antibody) at dilution of 1:500