

ZNF483 antibody



Catalog Number: 117182

Product name

ZNF483 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF483 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of ZNF483 recombinant protein (Accession Number: XM_047422863). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117182(ZNF483 Antibody) at dilution of 1:300