

ZZZ3 antibody



Catalog Number: 117275

Product name

ZZZ3 antibody

Specificity

Human; other species not tested.

Antibody description

ZZZ3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 150-160 kDa

Preparation

This antibody was obtained by immunization of ZZZ3 recombinant protein (Accession Number: XM_005270729). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117275(ZZZ3 antibody) at dilution of 1:1000