

DEF6 Antibody (C-term)



Catalog Number: 134868

General Information

Product name:

DEF6 Antibody (C-term)

Immunogen:

Applications:

WB, ELISA

Specificity:

Human

Antibody description:

Peptide Affinity Purified Rabbit Polyclonal Antibody (Pab)

Properties

Preparation:

This DEF6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 500-528 amino acids from the C-terminal region of human DEF6.

Formulation:

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is

purified through a protein A column, followed by peptide affinity purification.

Storage:

Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles. Store undiluted.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB: 1:1000

Validations



DEF6 Antibody (C-term) (Cat. #134868) western blot analysis in A375 cell line lysates (35ug/lane). This demonstrates the DEF6 antibody detected the DEF6 protein (arrow).