

Anti-EVI2B antibody



Catalog Number: 141082

Product name

Anti-EVI2B antibody

Specificity

Human EVI2B / CD361

Antibody description

Rabbit Polyclonal to Human EVI2B

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human EVI2B / CD361, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 5-20 µg/ml


IP: 1-4 µl/mg of lysate

Validations



Lanes	A	B
Items	Mouse heart tissue lysate	Mouse spleen tissue lysate
Sample Volume (µg/lane)	30	30
Gel	13% SDS-PAGE reducing gel	
Recommended Concentration	5-20 µg/ml	
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG [1:1], at 1:5000 dilution. Developed using Odyssey imaging system.	
Explanation	Predicted band size : 49 kDa Observed band size : 50 kDa	

Human EVI2B / CD361 Western blot (WB) 15574



Lanes	A	B
Items	HeLa	Jurkat
Sample (whole cell lysate)	0.5 mg	
IP antibody quantity	4 µl.	
Protein G agarose	15 µl of 50% Protein G Agarose	
Gel	13% SDS-PAGE reducing gel	
Primary antibody	KLH-SMCC- <i>hmoA2030</i> +KLH-SMCC- <i>hmoA2213</i> antibody at 10 µg/ml	
Secondary antibody	Clean-Blot5 IP Detection Reagent (HRP) at 1:500 dilution	

Human EVI2B / CD361 Immunoprecipitation(IP) 15575