

Anti-FCGR2A antibody



Catalog Number: 141206

Product name

Anti-FCGR2A antibody

Immunogen

[Human CD32a / Fc gamma RIIA / FCGR2A \(167 Arg, His Tag\) recombinant protein](#)

Specificity

Human CD32a/FCGR2A/Fc gamma RIIA

Antibody description

Mouse Monoclonal to Human FCGR2A

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD32a/FCGR2A/Fc gamma RIIA and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Human CD32a/FCGR2A/Fc gamma RIIA Flow Cytometry (FC) 15885