

# Anti-CD58 antibody



Catalog Number: 141242

## Product name

Anti-CD58 antibody

## Immunogen

[Human LFA-3 / CD58 \(His Tag\) recombinant protein](#)

## Specificity

Human LFA-3 / CD58

## Antibody description

Mouse Monoclonal to Human CD58

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human LFA-3 / CD58 and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



Human LFA-3 / CD58 Flow Cytometry (FC) 15900