

# Anti-FCGR2 antibody



Catalog Number: 141246

## Product name

Anti-FCGR2 antibody

## Immunogen

[Human CD32a / Fc gamma RIIA / FCGR2A \(167 Arg, His Tag\) recombinant protein](#)

## Specificity

Human CD32a

*Has cross-reactivity with Human CD32b*

## Antibody description

Mouse Monoclonal to Human FCGR2

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD32a/FCGR2A/Fc gamma RIIA and conjugated with PE under optimum conditions, the unreacted PE was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



Human CD32a/FCGR2A/Fc gamma RIIA Flow Cytometry (FC) 15890