

Anti-NRP2 antibody



Catalog Number: 141304

General Information

Product name:

Anti-NRP2 antibody

Immunogen:

[Human Neuropilin-2/NRP2 \(His Tag\) recombinant protein](#)

Applications:

FCM

Specificity:

Human NRP2 / Neuropilin-2

Antibody description:

Mouse Monoclonal to Human NRP2

Properties

Preparation:

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human NRP2 / Neuropilin-2. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation:

0.2 µm filtered solution in PBS

Storage:

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality:


Monoclonal

Ig Type: Mouse IgG1

Applications

FCM: 0.5-2 µg/test

Validations



Human NRP2/Neuropilin-2 Flow Cytometry (FC)
15128