

Anti-HNF4G antibody



Catalog Number: 141305

General Information

Product name:

Anti-HNF4G antibody

Immunogen:

Applications:

WB

Specificity:

Mouse HNF4G

Antibody description:

Rabbit Polyclonal to Mouse HNF4G

Properties

Preparation:

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Mouse HNF4G, and purified by antigen affinity chromatography.

Formulation:

0.2 µm filtered solution in PBS

Storage:

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB: 5-20 µg/ml

Validations



Items	Lanes	A
Sample		Mouse spleen tissue lysate
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-20 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution. Developed using Odyssey imaging system.
Explanation		Predicted band size : 46 kDa Observed band size : 46 kDa

HNF4G Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot