

# Anti-HNF4G antibody



Catalog Number: 141305

## Product name

Anti-HNF4G antibody

## Specificity

Mouse HNF4G

## Antibody description

Rabbit Polyclonal to Mouse HNF4G

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Mouse HNF4G, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-20 µg/ml

## Validations



Lanes	
Items	A
Sample	Mouse spleen tissue lysate
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 46 kDa Observed band size : 46 kDa

HNF4G Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot