

Anti-IGFBP6 antibody



Catalog Number: 141424

Product name

Anti-IGFBP6 antibody

Immunogen

[Mouse IGFBP6 / IBP6 \(His Tag\) recombinant protein](#)

Specificity

Mouse IGFBP6 / IBP6

Antibody description

Rabbit Monoclonal to Mouse IGFBP6

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IGFBP6 / IBP6 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse IGFBP6 / IBP6.

IHC-P: 1-10 µg/mL

Validations



Mouse IGFBP6 / IBP6 Immunohistochemistry(IHC)
15170