

Anti-FUT4 antibody



Catalog Number: 141455

Product name

Anti-FUT4 antibody

Specificity

Mouse FUT4

Antibody description

Rabbit Polyclonal to Mouse FUT4

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Mouse FUT4, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



| Items | Lanes | A |
|----------------------------|-------|--|
| Sample (whole cell lysate) | | HEK293 |
| Sample Volume (µg/lane) | | 30 |
| Gel | | 13% SDS-PAGE reducing gel |
| Recommended Concentration | | 5-20 µg/ml |
| Secondary Antibody | | DyLight 800-labeled Antibody to Rabbit IgG [H+], at 1:5000 dilution. |
| Explanation | | Developed using Odyssey imaging system. Predicted band size : 59 kDa Observed band size : 70 kDa |

Mouse FUT4 Western blot (WB) 14878