

Anti-Hdac7 antibody



Catalog Number: 141530

Product name

Anti-Hdac7 antibody

Specificity

Human Hdac7

Antibody description

Rabbit Polyclonal to Human Hdac7

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human Hdac7, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 4-8 µl/mg of lysate

Validations



Item	Lanes	A
Sample (whole cell lysate)		Hdca
Sample quantity		0.5 mg
IP antibody quantity		4 µL
Protein G agarose		15 µl of 50% Protein G agarose
Gel		13% SDS-PAGE reducing gel
Primary antibody		KLH-SMCC-SINOA734+KLH-SMCC-SINOA769 antibody at 10 µg/ml
Secondary antibody		Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution

Hdac7 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation