

# Anti-Hdac7 antibody



Catalog Number: 141530

## Product name

Anti-Hdac7 antibody

## Specificity

Human Hdac7

## Antibody description

Rabbit Polyclonal to Human Hdac7

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human Hdac7, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IP

## Dilutions

IP: 4-8 µl/mg of lysate

## Validations



Hdac7 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation