

Anti-HA antibody



Catalog Number: 141570

Product name

Anti-HA antibody

Immunogen

[H9N2 HA \(His Tag\) recombinant protein](#)

Specificity

Influenza H9N2 (A/Chicken/Hong Kong/G9/97)
Hemagglutinin/ HA1

Antibody description

Rabbit Polyclonal to H9N2 HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza H9N2 (A/Chicken/Hong Kong/G9/97) Hemagglutinin/ HA1 . Influenza H9N2 (A/Chicken/Hong Kong/G9/97) Hemagglutinin/ HA1 specific IgG was purified by Influenza H9N2 (A/Chicken/Hong Kong/G9/97) Hemagglutinin/ HA1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, FCM, ICC/IF, IF, IP

Dilutions

WB: 0.5-5 µg/ml

Validations



Items	Lanes
Sample (Recombinant protein)	A H9N2-CH-97-HA1-His (Cat#40036-V08H1)
Sample Volume (ng/lane)	50
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	0.5-5 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution. <small>Developed using Odyssey imaging system.</small>

Influenza H9N2 (A/Chicken/Hong Kong/G9/97) Hemagglutinin/ HA1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot