

Anti-GFRAL antibody



Catalog Number: 141587

Product name

Anti-GFRAL antibody

Specificity

Human GFRAL

Antibody description

Rabbit Polyclonal to Human GFRAL

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GFRAL, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Items	Lanes
Sample (whole cell lysate)	A HeLa
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	DyLight 800-labeled Antibody to Rabbit IgG [H+], at 1:5000 dilution.
Explanation	Developed using Odyssey imaging system. Predicted band size : 45 kDa Observed band size : 45 kDa

Human GFRAL Western blot (WB) 15450