

Anti-PROCR antibody



Catalog Number: 141704

Product name

Anti-PROCR antibody

Immunogen

[Mouse Epcr/PROCR \(His Tag\) recombinant protein](#)

Specificity

Mouse Epcr/PROCR

Antibody description

Rabbit Monoclonal to Mouse PROCR

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Epcr/PROCR and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



Mouse Epcr/PROCR Flow Cytometry (FC) 15367