

# Anti-HLA-DR antibody

Catalog Number: 141778

## Product name

Anti-HLA-DR antibody

## Specificity

Human HLA-DR

## Antibody description

Rabbit Polyclonal to Human HLA-DR

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human HLA-DR, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P


## Dilutions

IHC-P: 0.1-2µg/mL

## Validations



Human HLA-DR Immunohistochemistry(IHC)



15451

Human HLA-DR Immunohistochemistry(IHC)

15452