

Anti-FOXK1 antibody



Catalog Number: 141797

Product name

Anti-FOXK1 antibody

Specificity

Human FOXK1

Antibody description

Rabbit Polyclonal to Human FOXK1

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FOXK1, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC-P, ICC/IF, IF

Dilutions


IHC-P: 0.1-2 µg/mL

ICC/IF: 0.5-3µg/mL


Validations



Human FOXK1 Immunofluorescence(IF) 15393



Human FOXK1 Immunohistochemistry(IHC)



15392

Human FOXK1 Immunohistochemistry(IHC)
15391