

Anti-TET1 antibody



Catalog Number: 141815

Product name

Anti-TET1 antibody

Specificity

Human TET1

Antibody description

Rabbit Polyclonal to Human TET1

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human TET1, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 1-5 μ g/mL

Validations



Human TET1 Immunohistochemistry(IHC) 15447