

Anti-DSC2 antibody



Catalog Number: 141892

Product name

Anti-DSC2 antibody

Immunogen

[Rat DSC2 \(His Tag\) recombinant protein](#)

Specificity

Rat DSC2

Antibody description

Rabbit Monoclonal to Rat DSC2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat DSC2 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.5-5 µg/mL

Validations



Rat DSC2 Immunohistochemistry(IHC) 14556