

Anti-ING1 antibody



Catalog Number: 141970

Product name

Anti-ING1 antibody

Specificity

Mouse ING1

Antibody description

Rabbit Polyclonal to Mouse ING1

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Mouse ING1, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Lanes	A	B	C	D
Items				
Sample (whole cell lysate)	HepG2	HeLa	MCF7	Jurkat
Sample Volume (µg/lane)	30	30	30	30
Gel	13% SDS-PAGE reducing gel			
Recommended Concentration	5-20 µg/ml			
Secondary Antibody	DyLight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.			
	Developed using Odyssey imaging system.			
Explanation	Predicted band size : 47 kDa Observed band size : 55 kDa Additional bands at : 100 kDa and 33 kDa (We are unsure as to the identity of these extra bands.)			

ING1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot