

Anti-INS antibody



Catalog Number: 141977

Product name

Anti-INS antibody

Immunogen

[Human Insulin / INS recombinant protein](#)

Specificity

Human Insulin / INS

Antibody description

Rabbit Polyclonal to Human INS

Preparation

Produced in rabbits immunized with purified, recombinant Human Insulin / INS . Insulin / INS specific IgG was purified by Human Insulin / INS affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 µg/mL

Validations



Human Insulin / INS Immunohistochemistry(IHC)
15041