

# Anti-YWHAB antibody



Catalog Number: 142023

## Product name

Anti-YWHAB antibody

## Immunogen

[Cynomolgus 14-3-3 beta/YWHAB recombinant protein](#)

## Specificity

Cynomolgus 14-3-3 beta / YWHAB

## Antibody description

Rabbit Polyclonal to Cynomolgus YWHAB

## Preparation

Produced in rabbits immunized with purified, recombinant Cynomolgus 14-3-3 beta / YWHAB . 14-3-3 beta / YWHAB specific IgG was purified by Cynomolgus 14-3-3 beta / YWHAB affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IP

## Dilutions

IP: 1-4 µl/mg of lysate

## Validations



Item	Lanes	A	B
Sample (whole cell lysate)		HeLa	HepG2
Sample quantity		0.5 mg	
IP antibody quantity		4 µL	
Protein G agarose		15 µl of 50% Protein G Agarose	
Gel		15% SDS-PAGE reducing gel	
Primary antibody		Cyno YWHAB antibody at 10 µg/ml	
Secondary antibody		Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution	

Cynomolgus 14-3-3 beta/YWHAB Immunoprecipitation(IP) 14755