

Anti-Rock1 antibody



Catalog Number: 142072

Product name

Anti-Rock1 antibody

Specificity

Human Rock1

Antibody description

Rabbit Polyclonal to Human Rock1

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human Rock1, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 4-8 µl/mg of lysate

Validations



Items	A
Sample (whole cell lysate)	Jurkat
Sample quantity	0.5 mg
IP antibody quantity	4 µg
Protein G agarose	15 µl of 50% Protein G Agarose
Gel	15% SDS-PAGE reducing gel
Primary antibody	KLJ1-SMCC-sinoA2148+KLJ1-SMCC-sinoA2151 antibody at 10 µg/ml
Secondary antibody	Dylight 800-labeled antibody to Rabbit IgG (H+L) ₂ at 1:5000 dilution

Rock1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation