

Anti-DHRS9 antibody



Catalog Number: 142235

Product name

Anti-DHRS9 antibody

Immunogen

[Human DHRS9/3-alpha-HSD \(His Tag\) recombinant protein](#)

Specificity

Human DHRS9 / 3-alpha-HSD

Antibody description

Rabbit Polyclonal to Human DHRS9

Preparation

Produced in rabbits immunized with purified, recombinant Human DHRS9 / 3-alpha-HSD . DHRS9 / 3-alpha-HSD specific IgG was purified by Human DHRS9 / 3-alpha-HSD affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, IP

Dilutions


WB: 10-20 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human DHRS9 / 3-alpha-HSD.


IP: 1-4 uL/mg of lysate

Validations



Lanes	A
Items	
Sample (whole cell lysate)	THP1
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed	using Odyssey imaging system.
Explanation	Predicted band size : 35 kDa Observed band size : 35 kDa

DHRS9 / 3-alpha-HSD Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



Lanes	A
Items	
Sample (whole cell lysate)	THP1
Sample quantity	0.5 mg
IP antibody quantity	2 µg
Protein G agarose	15 µl of 50% Protein G Agarose
Gel	13% SDS-PAGE reducing gel
Primary antibody	His-DHRS9 antibody at 10 µg/ml
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution

DHRS9 / 3-alpha-HSD Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation